



Wipak[®] ECOKING

Soft and Noiseless Specialty Web



Wipak® ECOKING

Soft and Noiseless Specialty Web

Wipak® Ecoking is special polyethylene based film developed for urine and drainage applications. It is soft and noiseless and what is most important it is PVC-free.

Technical Details

- PE-based specialty web
- 80, 150 and 180 microns

Applications

- Suitable for urine and drainage bags
- For Irradiation sterilization

Benefits

- Comfortable in use due to soft and noiseless structure
- Environmental solution as no PVC is used
- Can be printed and heat sealed

For further Information please contact:
mdi@wipak.com
www.wipak.com



Involve. Inspire. Impress.

Wipak® is a registered trademark of Wipak Oy. Wipak reserves right to make changes at any time without notice.