



Wipak[®] NICE

Sustainable Replacement for Rigid Films



Wipak® NICE

Sustainable Replacement for Rigid Films

Wipak® NICE is semi rigid blown coextruded multilayer film intended in use as formable bottom web and replace rigid films. Suitable for EO and irradiation sterilization.

Technical Details

- Coextruded PE/PA/PE forming film
- 230 - 480 microns

Benefits

- Sustainable alternative with even 25% lighter packages
- Very stable, good formability and remaining corner thicknesses

Applications

- Sealable with coated and treated papers
- Ideal for OR kits, safety syringes and wide variety of special devices with sharp edges
- Replacing rigid films like APET, PETG and PVC
- For steam and EO sterilization

For further Information please contact:

mdi@wipak.com

www.wipak.com



Involve. Inspire. Impress.

Wipak® is a registered trademark of Wipak Oy. Wipak reserves right to make changes at any time without notice.