

# WINDOWBAG G2

Wipak's patented sterile barrier system concept



- Excellent protection of sterile medical devices and/or surgical draping packs
- Most convenient aseptic opening and sterile presentation
- Superior breathability for Ethylene Oxide sterilization



# WINDOWBAG G2

## Most convenient aseptic opening and sterile presentation

- WindowBag G2 gives wide peelable film-to-film window for opening and remaining the content
- No contamination risk from loose particles of coating adhesive
- The strong peelable seal facilitate clean, fibre-free, peel opening allowing aseptic presentation of the packed product

## Easy converting

- WindowBags can be easily closed with the impulse or rotary sealer
- The non peelable top of the bag prevents any seal damage during packing and transportation

## Superior breathability

- Use of uncoated Tyvek® provides superior breathability and does not limit the air and gas penetration / removal by overall coating
- Reduced risk of broken seals during sterilization process due to the position of the non peelable function of Tyvek®

## Excellent protection

- Window Bag G2 is designed to lower risk of broken seals
- Unique multilayer structure film is used with strong and flexible features
- Non peelable top seal avoids damages during transport and sterilization



Involve. Inspire. Impress.